

Cellulose-tethered tetradentate chelating agent for pyridoxine determination

Vera M. Ostrovskaya, Dmitry Yu. Marchenko and Alexander V. Ivanov

All standard solutions were prepared from chemically pure grade reagents in doubly distilled water. Pyridoxine hydrochloride [3,5-bis(hydroxymethyl)-2-methylpyridin-3-ol hydrochloride, Vitamine B₆, min 98%]; aniline (synthesis grade, 99%), phenol (synthesis grade, 99%, Merck) were used. Chromatographic paper containing 99% of hydrolyzed α -cellulose fibers (beats weight 80 g m⁻², thickness 0.18-0.20 mm, the content of Cu, Zn and Fe less than 0.5 ppm) was used as a carrier [S1].

The solid-phase reagent based on 2-hydroxy-4-nitroso-4-{2-[5-(2-diazo-5-nitrophenyl)-1,5-dioxapentyl]-4-nitrophenyl}-4-azabutylcellulose 4-toluenesulfonate (RIP-Diazo) in the form indicator strips was synthesized, prepared and azo coupled with pyridoxine to form of RIP-Azo-Pyridoxine as described [S2] (see Scheme 1 of the main text). Epoxy cellulose **1** was processed with 1,5-bis(2-amino-5-nitrophenyl)-1,5-dioxapentane to form of 2-hydroxy-3-{2-[5-(2-amino-5-nitrophenyl)-1,5-dioxapentyl]-4-nitrophenylamino}propylcellulose **2**, which was subjected to diazotization to form a reagent indicator paper RIP-Diazo **3**. Under the action of RIP-Diazo onto pyridoxine, a red-orange tetradentate reagent RIP-Azo-Pyridoxine **4** was formed. Test strips (TS) with a polymer holder 9.8 ± 0.05 mm and paper TS with a size of $40 \times 9.8 \pm 0.05$ mm of yellow-white color (TS RIP-Diazo) were prepared from RIP-Diazo for quantitative sorption-spectrophotometric determination of pyridoxine. The diffuse reflectance spectra of the test strips [in the 400-700 nm wavelength range, as the Kubelka-Munk function $F = (1-R)^2/2R$, where R is the reflection coefficient] were recorded on Eye-One Pro mini-spectrophotometer (X-Rite Inc., USA). The determination was based on azo coupling reaction with the TS RIP-Diazo to form test strips with azo-pyridoxine (TS RIP-Azo) and the subsequent complexation reaction of the tetradentate TS RIP-Azo with Cu^{II} ion, selected after testing for chromogenic reactions of pyridoxine with 14 metal ions: Ag, Al, Au^{III}, Ba, Bi^{III}, Fe^{III}, Ca, Cd, Co^{II}, Cu^{II}, Mg, Mn^{II}, Ni, Pb and Zn. To do this, aqueous solution of pyridoxine (0.5 ml) was adjusted to optimal pH 3–8 with a buffer solution. The TS RIP-Diazo with the holder was immersed with a paper end into the manufactured pyridoxine sample for 1–5 min, the reaction solution was removed, the standard metal solution (50 μ l) and an acetate buffer solution up to pH 5 were added. Color change of the TS was observed before and after the reaction with the metal. The diffuse reflectance spectrum of the dried test-strip RIP-Azo and the product of its reaction with the metal ion (TS RIP-M) were

recorded. To determine pyridoxine in a pharmaceutical formulation, 1 tablet (5 mg) or averaged sample weighing 5.00 mg was dissolved in water in a volumetric flask, diluted up to 100 ml, and 20 ml of sample was taken, placed in a Petri dish, the solution was adjusted up to pH of 5–7 and standard solution of Cu^{II} (50 μl) was added. The paper TS RIP-Diazo was lowered into the solution for 20 min, dried, and the diffuse reflectance spectrum of TS was registered and the function F was calculated. The concentration of the main substance in the sample was determined from the calibration curve for the range of pyridoxine concentrations of 0.2–10 mg dm^{-3} .

References

- [S1] M. I. Berezhnaya, N. P. Ivanova, V. M. Ostrovskaya, P. P. Gura, Yu. N. Prikhod'ko and P.V. Poddubny, *SU Patent 1219698*, 1986.
- [S2] V. M. Ostrovskaya, L. K. Shpigun and Ya. V. Shushenachev, *RU Patent 2680391*, 2019.